



Endobon® Xenograft Granules **A RegenerOss™ Product**

A Scaffold For Building Bone

- **Fully Deproteinated By Two-Step, High Temperature Manufacturing Process For Safety From Bacteria, Viruses And Prions**
- **Osseoconductive Properties Enable Bone To Grow Directly On The Ceramic Surface And Through The Entire Graft**
- **Excellent Handling Characteristics For Easy Transfer To The Defect Site**




Photography Courtesy Of: Dr. John Lupovici

The Elements For Building Bone

Comparable To Natural Bone:

- 1 To 3µm Crystallite Size
- HA Content >95% After Processing
- Interconnecting Micro And Macro Pores For Quick Bony Integration, Which May Lead To Primary Mechanical Stability And Vascular Ingrowth



Endobon® Xenograft Granules At 20x

Indications For Dental And/Or Oral Surgical Procedures:

- Aesthetic Contouring Defects
- Peri-Implant Defects
- Ridge Augmentation
- As A Graft Extender



Clinical Photograph Courtesy Of:
Dr. José Luis Calvo Guirado

Ordering Information:

Endobon Xenograft Granules - A RegenerOss™ Product

Catalog Number	Volume	Particle Size
ROX05	0.5cc	500-1000µm
ROX10	1.0cc	500-1000µm
ROX20	2.0cc	500-1000µm

Distributed By:



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside The U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com

Call Your BIOMET 3i Sales Representative Today

Sign Up For BIOMET 3i's Electronic Newsletter "BIOMET 3innovations."
Simply Go Online To www.biomet3i.com

Responsible Manufacturer: BIOMET Orthopaedics Switzerland GmbH

RegenerOss is a trademark of BIOMET 3i, Inc. BIOMET is a registered trademark and BIOMET 3i and design are trademarks of BIOMET, Inc.
©2008 BIOMET 3i, Inc. All rights reserved.

Endobon is a registered trademark of BIOMET Deutschland GmbH.



ART1030

REV A 01/08

SUBSIDIARIES

AUSTRALIA
Phone: +61-2-9855-4444
Fax: +61-2-9888-9900

AUSTRIA
Phone: +43-(0)6235-200-45
Fax: +43-(0)6235-200-45-9

BELGIUM
Phone: +32-2-5410290
Fax: +32-2-5410291

BRAZIL
Phone: +55-11-5081-4405
Fax: +55-11-5081-7484

CANADA
Phone: +514-956-9843
Fax: +514-956-9844

FRANCE
Phone: +33-1-41054343
Fax: +33-1-41054340

GERMANY
Phone: +49-721-255177-10
Fax: +49-721-255177-73

IRELAND
Phone: +35-31-477-3925
Fax: +35-31-402-9590

MEXICO
Phone: +52-55-5679-1619
Fax: +52-55-5684-8098

THE NETHERLANDS
Phone: +31-(0)78-629-2800
Fax: +31-(0)78-629-2801

NEW ZEALAND
Phone: +64-508-122-221
Fax: +64-508-133-331

NORDIC REGION
Phone: +46-40-17-6090
Fax: +46-40-17-6099

PORTUGAL
Phone: +351-21-000-1645
Fax: +351-21-000-1675

SPAIN
Phone: +34-93-470-59-50
Fax: +34-93-372-11-25

SWITZERLAND
Phone: +41-44-380-46-46
Fax: +41-44-383-46 55

U.K.
Phone: +44-1628-829314
Fax: +44-1628-820182

DISTRIBUTORS

ARGENTINA
Dentalmax, SA
Phone: +541-1482-71001
Fax: +541-1482-67373

CHILE
Cybel, SA
Phone: +56-2-2321883
Fax: +56-2-2330176

CHINA
Atek Inc.
Phone: +86-21-6329-1265
Fax: +86-21-6329-1620

COLOMBIA
3i Colombia
Phone: +571-612-9362
Fax: +571-620-6412

COSTA RICA
Implantec S.A.
Phone: +506-2-256411
Fax: +506-2-247620

EL SALVADOR
Dentimer SA de CV
Phone: +503-263-6350
Fax: +503-263-6676

GREECE
Kostas Kornisortis and Co.
Phone: +30-2310-501-651
Fax: +30-2310-522-417

HONG KONG
Osttek Inc., Ltd.
Phone: +852-8121-6601
Fax: +852-3747-3754

ISRAEL
H.A. Systems
Phone: +972-3-6138777
Fax: +972-3-6138778

ITALY
Biomax, srl.
Phone: +39-0444-913410
Fax: +39-0444-913695

JAPAN
Implant Innovations Japan
Phone: +81-66-868-3012
Fax: +81-66-868-2444

KOREA
Jungsan Biomed Corp.
Phone: +82-2-516-1808
Fax: +82-2-514-9434

LEBANON
3i MENA s.a.l.
Middle East And North Africa
Phone: +961-1-694000
Fax: +961-1-694222

PARAGUAY
Andres H. Arce y Cia SRL
Phone: +595-21-208185
Fax: +595-21-496291

POLAND
Dental Depot
Phone: +48-71-341-3091
Fax: +48-71-343-6560

RUSSIA
Com-Dental
Phone: +7-495-797-6686
Fax: +7-499-242-9567

SINGAPORE
Asia Implant Support & Services
Phone: +65-6223-2229
Fax: +65-6220-3538

SOUTH AFRICA
Selective Surgical CC
Phone: +27-11-991-7007
Fax: +27-11-672-1391

TAIWAN
Kuo Hwa Dental Suppliers Co., Ltd.
Phone: +886-2-2226-1770
Fax: +886-2-2226-8747

THAILAND
3i (Thailand) Co., LTD.
Phone: +662-252-6685
Fax: +662-252-6686

UKRAINE
Com-Dental
Phone: +38-067-7007667
Fax: +38-044-5017117

URUGUAY
Pro3implant S.R.L.
Phone: +598-2-4034163
Fax: +598-2-4034163