

**Straumann® BLX Implant System System Overview**

**BLX PROSTHETIC**

**Impression Posts  
Implant Analogs  
Digital Impression Posts  
Repositionable Analogs**

**Temporary Abutments**



AH = abutment height,  
GH = gingiva height,  
Ø = diameter


**Variobase®**



AH = abutment height,  
GH = gingiva height,  
Ø = diameter

**Hybrid Abutments**


**Titanium Abutments**



AH = abutment height,  
GH = gingiva height,  
Ø = diameter

**RB/WB prosthetics – fit on RB and WB implants (Ø 3.5 mm - Ø 6.5 mm)**

- crowns with slim emergence profiles
- bridges




**Temporary Abutments for Crown, AH 11 mm**

GH	Ø 3.8 mm (RB/WB) Ø 4.5 mm (RB/WB) Ø 6.0 mm (RB/WB)		
	1.5 mm	064.4362	064.4372
2.5 mm	064.4363	064.4373	064.4382
3.5 mm	064.4364	064.4374	064.4383

**Straumann® BLX Implant TorcFit™ Connection**

**WB prosthetics – fit only on WB implants (Ø 5.0 mm - Ø 6.5 mm)**

- crowns with wide emergence profiles



**Temporary Abutments for Crown, AH 11 mm**

GH	Ø 5.5 mm (WB)		Basal Screw
	0.75 mm	064.4391	
1.5 mm	064.4392		
2.5 mm	–		
3.5 mm	–		

**Temporary Abutments for Bridge, AH 5.0 mm**

GH	Ø 4.5 mm (RB/WB)
1.5 mm	064.4352

**VB Basal Screw**

065.0036
----------

**Variobase® for Crown, AH 5.5 mm**

GH	Ø 3.8 mm (RB/WB) Ø 4.5 mm (RB/WB)	
	1.5 mm	062.4934
2.5 mm	062.4935	062.4945
3.5 mm	062.4936	062.4946

**VB Basal Screw**

065.0036
----------

**Variobase® for Crown AS, AH 5.5 mm**

GH	Ø 4.5 mm (RB/WB)
1.5 mm	062.4972

**Variobase® C for Sirona® Cerec**

Ø	Ø 3.8 mm Ø 4.5 mm
3.8 mm	062.4981*
4.5 mm	062.4982*

**Basal Screw VB AS**

065.0037
----------

**Variobase® for Bridge/Bar Cylindrical, AH 3.5 mm**

GH	Ø 4.5 mm (RB/WB)
1.5 mm	062.4961

**Burn-out Coping for Variobase® for Bridge/Bar Cylindrical**

Ø 4.5 mm (RB/WB)
065.0017

**VB Basal Screw**

065.0036
----------

**Variobase® for Crown, AH 5.5 mm**

GH	Ø 5.5 mm (WB)	
	0.75 mm	062.4953
1.5 mm	062.4954	
2.5 mm	–	
3.5 mm	–	

**VB Basal Screw**

065.0036
----------

**Variobase® for Crown AS, AH 5.5 mm**

GH	Ø 5.5 mm (WB)
1.5 mm	062.4971

**VB Basal Screw AS**

065.0037
----------

**Burn-out Coping for Variobase® for Crown**

Ø 5.5 mm (WB)
065.0016

**Burn-out Coping for Variobase® for Crown AS**

Ø 5.5 mm (WB)
065.0019

**SRA Screw-retained Abutments**

Angle	Ø 4.6 mm (RB/WB)		
	0°	–	–
17°	–	–	–
30°	–	–	–

**SRA Closed-tray Impression for Screw-retained Abutments**

Ø 4.6 mm	Crown	
	025.2246	025.2244

**SRA Open-tray Impression for Screw-retained Abutments**

Ø 4.6 mm	Crown	
	025.2244	025.0012

**SRA Burn-out Copings, Plastic**

Ø 4.6 mm	Crown	
	023.4748	023.4758

**SRA Titanium Copings Regular**

Ø 4.6 mm	Crown	
	023.4747	023.4751

**SRA Titanium Copings Slim**

Ø 4.6 mm	Crown	
	024.0023	024.0024

**SRA Variobase® for Bridge/Bar Cylindrical**

Ø 4.6 mm	Crown	
	023.0028	023.4763

**SRA Burn-out Coping for Variobase® for Bridge/Bar Cylindrical**

Ø 4.6 mm	023.0032
----------	----------

**SRA Occlusal Screw**

023.4763
----------

**SRA Lab Processing Screw**

L 20mm	025.0006
L 10mm	025.0052

**RB/WB Pre-milled Abutment Blank**

Ø	11.5mm	Medentika Holder	062.4601
12mm	Amann Girrbach	062.4603	

**Novaloc®**

Angle	Ø 3.8 mm (RB/WB)	
	0°	15°
1.5 mm	062.4501	–
2.5 mm	062.4502	062.4507
3.5 mm	062.4503	062.4508
4.5 mm	062.4504	062.4509
5.5 mm	062.4505	062.4510
6.5 mm	062.4506	062.4511
7.5 mm	–	062.4512

**Processing Package titanium**

2010.601-STM
--------------

**Processing Package PEEK**

2010.611-STM
--------------

**Abutment removal kit**

RB/WB Abutment-Removal Screw		for removal of damaged abutments, TAN	065.0007
Removal Tool		for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0008
		for BLX Basal Screw, left-hand, length 21 mm, stainless steel	065.0009

**מפתח פרוטזי**

046.400		15 mm
046.401		21 mm
046.402		27 mm
046.119		
046.049		