

**Straumann®
BLX Implant System
System Overview**

**BLX
SURGERY**

Roxolid® SLA®

Roxolid® SLActive®

Implant Drivers

**Closure Cap
H 0.4 mm**

Healing Abutment

**Auxiliary
Instruments**

∅ 3.5 mm

∅ 3.75 mm

∅ 4.0 mm

∅ 4.5 mm

∅ 5.5 mm

∅ 6.5 mm

BLX
RB (Regular Base)
∅ 3.5 mm



BLX
RB (Regular Base)
∅ 3.75 mm



BLX
RB (Regular Base)
∅ 4.0 mm



BLX
RB (Regular Base)
∅ 4.5 mm




BLX
WB (Wide Base)
∅ 5.5 mm



BLX
WB (Wide Base)
∅ 6.5 mm






Length	∅ 3.5 mm	∅ 3.75 mm	∅ 4.0 mm	∅ 4.5 mm	∅ 5.5 mm	∅ 6.5 mm
6 mm	-	6 mm 061.4506	6 mm 061.5506	6 mm 061.6506		
8 mm	061.3508	8 mm 061.4508	8 mm 061.5508	8 mm 061.6508		
10 mm	061.3510	10 mm 061.4510	10 mm 061.5510	10 mm 061.6510		
12 mm	061.3512	12 mm 061.4512	12 mm 061.5512	12 mm 061.6512		
14 mm	061.3514	14 mm 061.4514	14 mm 061.5514	14 mm 061.6514		
16 mm	061.3516	16 mm 061.4516	16 mm 061.5516	16 mm 061.6516		
18 mm	061.3518	18 mm 061.4518	18 mm 061.5518	18 mm 061.6518		
6 mm	-	6 mm 061.4306	6 mm 061.5306	6 mm 061.6306	6 mm 061.8306	6 mm 061.9306
8 mm	061.3308	8 mm 061.4308	8 mm 061.5308	8 mm 061.6308	8 mm 061.8308	8 mm 061.9308
10 mm	061.3310	10 mm 061.4310	10 mm 061.5310	10 mm 061.6310	10 mm 061.8310	10 mm 061.9310
12 mm	061.3312	12 mm 061.4312	12 mm 061.5312	12 mm 061.6312	12 mm 061.8312	12 mm 061.9312
14 mm	061.3314	14 mm 061.4314	14 mm 061.5314	14 mm 061.6314		
16 mm	061.3316	16 mm 061.4316	16 mm 061.5316	16 mm 061.6316		
18 mm	061.3318	18 mm 061.4318	18 mm 061.5318	18 mm 061.6318		



RB (Regular Base)

RB (Regular Base)

WB (Wide Base)



Healing Abutments for Crown

AH 2 mm	Diameter (∅)			
	GH	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm		064.4202S	064.4212S	064.4222S
2.5 mm		064.4204S	064.4214S	—
3.5 mm		064.4206S	064.4216S	—

Healing Abutments for Crown

AH 4 mm	Diameter (∅)			
	GH	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm		064.4203S	064.4213S	064.4223S

Healing Abutment for Bridge RB/WB

AH 2 mm	Diameter (∅)
	GH
	∅ 5.0 mm
1.5 mm	064.4232S

Healing Abutment for Bridge RB/WB

AH 4 mm	Diameter (∅)
	GH
	∅ 5.0 mm
1.5 mm	064.4233S

Healing Abutments for Crown Wide Profile

AH 2 mm	Diameter (∅)		
	GH	∅ 6.0 mm	∅ 7.5 mm
0.75 mm		064.8201S	064.8203S

Healing Abutments for Crown Wide Profile

AH 4 mm	Diameter (∅)		
	GH	∅ 6.0 mm	∅ 7.5 mm
0.75 mm		064.8202S	064.8204S


Torque Control Device		066.1100
Ratchet		046.119

Holding Key		046.064
Implant depth gauge		066.2000

BL Bone Profiler 1	∅ 5.2		026.0022
BL Bone Profiler 2	∅ 6.6		026.0023
BL Bone Profiler 3	∅ 6.0		026.0024
Guiding Pin for BLX Implants			066.0025S

Implant-abutment compatibility

RB/WB and WB		RB Implants ∅ 3.5 - 4.5 mm	WB Implants ∅ 5.0 - 6.5 mm
«RB/WB» RB/WB Variobase® for Crown	 compatible with		
«WB» WB Variobase® for Crown	 compatible with		



**No WB abutments
on RB implants!**

International Headquarters
Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
+41 (0)61 965 11 11
+41 (0)61 965 11 01
www.straumann.com

Phone
Fax

© Institut Straumann AG, 2019. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.