


Straumann® BLX Implant System System Overview


BLX PROSTHETIC

**Impression Posts
Implant Analogs
Digital Impression Posts
Repositionable Analogs**

Temporary Abutments




Variobase®




Hybrid Abutments

Titanium Abutments



RB/WB prosthetics – fit on RB and WB implants (Ø 3.5 mm - Ø 6.5 mm)

- crowns with slim emergence profiles
- bridges




Temporary Abutments for Crown, AH 11 mm

GH	Ø 3.8 mm (RB/WB) Ø 4.5 mm (RB/WB) Ø 6.0 mm (RB/WB)		
	1.5 mm	064.4362	064.4372
2.5 mm	064.4363	064.4373	064.4382
3.5 mm	064.4364	064.4374	064.4383

Straumann® BLX Implant TorcFit™ Connection

WB prosthetics – fit only on WB implants (Ø 5.0 mm - Ø 6.5 mm)

- crowns with wide emergence profiles



Temporary Abutments for Crown, AH 11 mm

GH	Ø 5.5 mm (WB)		Basal Screw
	0.75 mm	064.4391	
1.5 mm	064.4392		
2.5 mm	–		
3.5 mm	–		

Variobase® for Crown, AH 5.5 mm

GH	Ø 3.8 mm (RB/WB) Ø 4.5 mm (RB/WB)		
	1.5 mm	062.4934	062.4944
2.5 mm	062.4935	062.4945	
3.5 mm	062.4936	062.4946	

VB Basal Screw

065.0036

Variobase® for Crown AS, AH 5.5 mm

GH	Ø 4.5 mm (RB/WB)	
	1.5 mm	062.4972

Variobase® C for Sirona® Cerec

Ø	Ø 4.5 mm (RB/WB)	
	3.8 mm	062.4981*
4.5 mm	062.4982*	

Basal Screw VB AS

065.0037

Variobase® for Bridge/Bar Cylindrical, AH 3.5 mm

GH	Ø 4.5 mm (RB/WB)	
	1.5 mm	062.4961

Burn-out Coping for Variobase® for Bridge/Bar Cylindrical

Ø 4.5 mm (RB/WB)	
065.0017	

VB Basal Screw

065.0036

Variobase® for Crown, AH 5.5 mm

GH	Ø 5.5 mm (WB)	
	0.75 mm	062.4953
1.5 mm	062.4954	
2.5 mm	–	
3.5 mm	–	

VB Basal Screw

065.0036

Variobase® for Crown AS, AH 5.5 mm

GH	Ø 5.5 mm (WB)	
	1.5 mm	062.4971

VB Basal Screw AS

065.0037

Burn-out Coping for Variobase® for Crown

Ø 5.5 mm (WB)	
065.0016	

Burn-out Coping for Variobase® for Crown AS

Ø 5.5 mm (WB)	
065.0019	

SRA Screw-retained Abutments

Angle	Ø 4.6 mm (RB/WB)		
	0°	–	–
17°	–	–	–
30°	–	–	–

SRA Closed-tray Impression for Screw-retained Abutments

Ø	Ø 4.6 mm	
	Crown	025.2246
Bridge/Bar	025.0014	

SRA Open-tray Impression for Screw-retained Abutments

Ø	Ø 4.6 mm	
	Crown	025.2244
Bridge/Bar	025.0012	

SRA Burn-out Copings, Plastic

Ø	Ø 4.6 mm	
	Crown	023.4748
Bridge	023.4758	

SRA Titanium Copings Regular

Ø	Ø 4.6 mm	
	Crown	023.4747
Bridge	023.4751	

SRA Titanium Copings Slim

Ø	Ø 4.6 mm	
	Crown	024.0023
Bridge	024.0024	

SRA Variobase® for Bridge/Bar Cylindrical

Ø	Ø 4.6 mm	
	Crown	023.0028
Bridge	incl. screw 023.4763	

SRA Occlusal Screw

023.4763

SRA Lab Processing Screw

L 20mm	025.0006
L 10mm	025.0052

RB/WB Pre-milled Abutment Blank

Ø	Ø 4.6 mm	
	11.5mm	Medentika Holder
12mm	Amann Girrbach	062.4603

incl. Screw 065.0036

WB Variobase® C for Sirona®

Ø	Ø 5.5 mm (WB)	
	5.5mm	062.4983*

RB/WB Pre-milled Abutment Blank

Ø	Ø 4.6 mm	
	11.5mm	Medentika Holder
12mm	Amann Girrbach	062.4607

incl. Screw 065.0036

SRA Analogs

Abutment Analog	023.4756
Repositionable Analog	025.0008

Straumann® CARES® Scanbody for Screw-retained Abutments

mm 4.6 Ø	025.0001
----------	----------

Anatomic Abutments

Angle	Ø 6.0 mm (RB/WB)	
	0°	–
17°	–	

GH	Ø 6.0 mm (RB/WB)	
	2.5 mm	062.4103
3.5 mm	062.4104	062.4154

Anatomic Abutment Basal Screw	065.0036
-------------------------------	----------

Novaloc®

Angle	Ø 3.8 mm (RB/WB)	
	0°	–
15°	–	

GH	Ø 3.8 mm (RB/WB)	
	1.5 mm	062.4501
2.5 mm	062.4502	062.4507
3.5 mm	062.4503	062.4508
4.5 mm	062.4504	062.4509
5.5 mm	062.4505	062.4510
6.5 mm	062.4506	062.4511
7.5 mm	–	062.4512

Processing Package titanium

2010.601-STM

Processing Package PEEK

2010.611-STM

Abutment removal kit

RB/WB Abutment-Removal Screw	for removal of damaged abutments, TAN	065.0007
Removal Tool	for BLX Basal Screw, left-hand, length 27 mm, stainless steel	065.0008
	for BLX Basal Screw, left-hand, length 21 mm, stainless steel	065.0009

מפתח פרוטזי

046.400	15 mm
046.401	21 mm
046.402	27 mm
046.119	
046.049	

**Straumann®
BLX Implant System
System Overview**

**BLX
SURGERY**

Roxolid® SLA®

Roxolid® SLActive®

Implant Drivers

**Closure Cap
H 0.4 mm**

Healing Abutment

**Auxiliary
Instruments**

∅ 3.5 mm

∅ 3.75 mm

∅ 4.0 mm

∅ 4.5 mm

∅ 5.5 mm

∅ 6.5 mm

BLX
RB (Regular Base)
∅ 3.5 mm



BLX
RB (Regular Base)
∅ 3.75 mm



BLX
RB (Regular Base)
∅ 4.0 mm



BLX
RB (Regular Base)
∅ 4.5 mm








BLX
WB (Wide Base)
∅ 5.5 mm



BLX
WB (Wide Base)
∅ 6.5 mm






Length	∅ 3.5 mm	∅ 3.75 mm	∅ 4.0 mm	∅ 4.5 mm	∅ 5.5 mm	∅ 6.5 mm
6 mm	-	061.4506	061.5506	061.6506		
8 mm	061.3508	061.4508	061.5508	061.6508		
10 mm	061.3510	061.4510	061.5510	061.6510		
12 mm	061.3512	061.4512	061.5512	061.6512		
14 mm	061.3514	061.4514	061.5514	061.6514		
16 mm	061.3516	061.4516	061.5516	061.6516		
18 mm	061.3518	061.4518	061.5518	061.6518		
6 mm	-	061.4306	061.5306	061.6306	061.8306	061.9306
8 mm	061.3308	061.4308	061.5308	061.6308	061.8308	061.9308
10 mm	061.3310	061.4310	061.5310	061.6310	061.8310	061.9310
12 mm	061.3312	061.4312	061.5312	061.6312	061.8312	061.9312
14 mm	061.3314	061.4314	061.5314	061.6314		
16 mm	061.3316	061.4316	061.5316	061.6316		
18 mm	061.3318	061.4318	061.5318	061.6318		



RB (Regular Base)

RB (Regular Base)

WB (Wide Base)



Healing Abutments for Crown

AH 2 mm	Diameter (∅)			
	GH	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm		064.4202S	064.4212S	064.4222S
2.5 mm		064.4204S	064.4214S	—
3.5 mm		064.4206S	064.4216S	—

Healing Abutments for Crown

AH 4 mm	Diameter (∅)			
	GH	∅ 4.0 mm	∅ 5.0 mm	∅ 6.5 mm
1.5 mm		064.4203S	064.4213S	064.4223S

Healing Abutment for Bridge RB/WB

AH 2 mm	Diameter (∅)
	GH
1.5 mm	064.4232S

Healing Abutment for Bridge RB/WB

AH 4 mm	Diameter (∅)
	GH
1.5 mm	064.4233S

Healing Abutments for Crown Wide Profile

AH 2 mm	Diameter (∅)		
	GH	∅ 6.0 mm	∅ 7.5 mm
0.75 mm		064.8201S	064.8203S

Healing Abutments for Crown Wide Profile

AH 4 mm	Diameter (∅)		
	GH	∅ 6.0 mm	∅ 7.5 mm
0.75 mm		064.8202S	064.8204S


Torque Control Device		066.1100
Ratchet		046.119

Holding Key		046.064
Implant depth gauge		066.2000

BL Bone Profiler 1	∅ 5.2		026.0022
BL Bone Profiler 2	∅ 6.6		026.0023
BL Bone Profiler 3	∅ 6.0		026.0024
Guiding Pin for BLX Implants			066.0025S

Implant-abutment compatibility

RB/WB and WB		RB Implants ∅ 3.5 - 4.5 mm	WB Implants ∅ 5.0 - 6.5 mm
«RB/WB» RB/WB Variobase® for Crown	 compatible with		
«WB» WB Variobase® for Crown	 compatible with		



**No WB abutments
on RB implants!**

International Headquarters
Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
+41 (0)61 965 11 11
+41 (0)61 965 11 01
www.straumann.com

Phone
Fax

© Institut Straumann AG, 2019. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.