


Straumann® Dental Implant System System Overview TISSUE LEVEL SURGERY MRK.034.Rev.04	Ø 3.3 mm				Ø 4.1 mm				Ø 4.8 mm														
		Standard Plus (SP) RN Ø 4.8 mm  Ø 3.3 mm				Standard Plus (SP) RN Ø 4.8 mm  Ø 4.1 mm				Standard Plus (SP) RN Ø 4.8 mm  Ø 4.8 mm			Standard Plus (SP) WN Ø 6.5 mm  Ø 4.8 mm										
Titanium SLA®		4 mm	–	Roxolid® SLA® only		4 mm	–		4 mm	–		4 mm	–										
		6 mm	–			6 mm	043.9255		6 mm	043.9395		6 mm	043.9485										
		8 mm	043.2605			8 mm	043.9265		8 mm	043.9405		8 mm	043.9495										
		10 mm	043.2615			10 mm	043.9275		10 mm	043.9415		10 mm	043.9505										
		12 mm	043.2625			12 mm	043.9285		12 mm	043.9425		12 mm	043.9515										
		14 mm	043.2635			14 mm	043.9295		14 mm	043.9435		14 mm	–										
Roxolid® SLActive®		6 mm	–			4 mm	033.0435		4 mm	033.0445		4 mm	033.0455										
		8 mm	033.5115			6 mm	033.5605		6 mm	033.5905		6 mm	033.6105										
		10 mm	033.5125			8 mm	033.5615		8 mm	033.5915		8 mm	033.6115										
		12 mm	033.5135			10 mm	033.5625		10 mm	033.5925		10 mm	033.6125										
		14 mm	033.5145			12 mm	033.5635		12 mm	033.5935		12 mm	033.6135										
Round Burs	 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004				 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004				 Round Bur Ø 2.3 mm 044.003 →  Round Bur Ø 3.1 mm 044.004														
Pilot and Twist Drills	 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill Pro Ø 2.8 mm 044.787 short 044.789 long				 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill Ø 2.8 mm 044.787 short 044.789 long →  Twist Drill PRO Ø 3.5 mm 044.791 short 044.793 long				 Pilot Drill Ø 2.2 mm 044.783 short 044.785 long →  Twist Drill PRO Ø 2.8 mm 044.787 short 044.789 long →  Twist Drill PRO Ø 3.5 mm 044.791 short 044.793 long →  Twist Drill PRO Ø 4.2 mm 044.795 short 044.797 long														
Alignment Pins Depth Gauges	 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 22.8 mm 046.800				 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 2.8 mm 046.800 →  Depth Gauge Ø 3.5 mm 046.802				 Alignment Pin Ø 2.2 mm 046.799 →  Depth Gauge Ø 22.8 mm 046.800 →  Depth Gauge Ø 3.5 mm 046.802 →  Depth Gauge Ø 4.2 mm 046.804														
Profile Drills	 SP Profile Drill 034.332				 SP Profile Drill 034.333				 SP Profile Drill 044.084														
Taps	 SP Tap Ø 3.3 mm 034.351				 SP Tap Ø 4.1 mm 034.352  BL Tap Ø 4.1 mm 034.349 רק לשתל באורך 4 ס"מ				 SP Tap Ø 4.8 mm 034.353  BL Tap Ø 4.1 mm 034.350 רק לשתל באורך 4 ס"מ														
Closure Caps Healing Caps	RN																						
	Closure Cap - Flat								Healing Abutments														
 048.3715  Ø 6 mm H 1.5 mm 048.3735  Ø 5.5 mm H 2 mm 048.0335  Ø 5.5 mm H 3 mm 048.0345  Ø 5.5 mm H 4.5 mm 048.0375																 H 0 mm 048.3755  H 2 mm 048.0385  H 3 mm 048.0395  H 4.5 mm 048.0535							