

Slim emergence profile (RB/WB) prosthetics

fit on all BLC and BLX implants (Ø 3.3 mm - Ø 6.5 mm)

- crowns with slim emergence profiles
- bridges

Implant-abutment compatibility

		RB Implants	WB Implants
«RB/WB» RB/WB Variobase® XC for Crown			
«WB» WB Variobase® XC for Crown			

Wide Emergence Profile (WB) prosthetics

fit only on WB implants (BLC Ø 4.5 mm - Ø 6.5 mm, BLX Ø 5.0 mm - Ø 6.5 mm)

- crowns with wide emergence profiles

Abutment Level, on Screw-retained Abutments, Ø 4.6 mm

Implant Analog for Screw-retained Abutments, Ø 4.6 mm, 0°	Open-tray Impression for Screw-retained Abutments, Ø 4.6 mm	Closed-tray Impression for Screw-retained Abutments, Ø 4.6 mm
025.0050*	Crown 025.2244 Bridge/Bar 025.0012	Crown 025.2246 Bridge/Bar 025.0014

Scanbody for Screw-retained Abutments, Ø 4.6 mm

Repositionable Implant Analog for Screw-retained Abutments, Ø 4.6 mm

Scanbody / RevEX SRA 4.6

025.0081

025.0008

025.0083

Implant Level

Analog	Open-tray	Closed-tray
065.0021	065.0031	065.0033
		065.4310

Straumann® Scanbody

Repositionable Implant Analog

065.0103

065.0023

Temporary Titanium Coping for Screw-retained Abutments, Ø 4.6 mm

Titanium Copings	Ø 4.6 mm (RB/WB)
	Crown 023.4747
	Bridge 023.4751
incl. screw	023.4763

Repositionable Implant Analog

023.4763

RB/WB Temporary Abutment, VITA CAD-Temp®

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
GH 1.5 mm	064.4361
	064.4371

Temporary Abutments for Crown, AH 10 mm

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)	Ø 6.0 mm (RB/WB)
GH 1.5 mm	064.4362	064.4372
	064.4363	064.4382
GH 2.5 mm	064.4374	064.4383
GH 3.5 mm	064.4364	064.4383
	065.0036	

Replacement screw

065.0036

Temporary Abutments for Bridge, AH 10 mm

Ø 4.5 mm (RB/WB)	
GH 1.5 mm	064.4352

Replacement screw

065.0036

Temporary Abutments for Crown, AH 11 mm

Ø 5.5 mm (WB)	
GH 0.75 mm	064.4391
	064.4392
GH 1.5 mm	064.4390

Replacement screw

065.0036

WB Temporary Abutment, VITA CAD-Temp®

Ø 5.5 mm (WB)	
GH 1.5 mm	064.4390

Screw-retained Abutments

Angle	0°	17°	30°
0.75 mm	—	—	—
1.5 mm	062.47225	—	—
2.5 mm	062.47235	—	—
3.5 mm	062.47245	062.47335	062.47435
4.5 mm	062.47255	062.47345	062.47445
5.5 mm	—	062.47355	062.47455

Replacement screw

065.0036

Variobase® for Bridge/Bar

Ø 4.6 mm
Crown 023.4748
Bridge 023.4758
incl. screw 023.4763

Lab Processing Screw

L 20mm	025.0006
L 10mm	025.0052

RB/WB SRA PLAN Abutment, POM

0°	17°	30°
025.0073V4	025.0074V4	025.0075V4

Occlusal Screw

023.4763

Occlusal Screw AS

025.0054

Lab Processing Screw

L 20mm	025.0006
L 10mm <td>025.0052</td>	025.0052

RB/WB Pre-milled Abutment Blanks

For MEDENTIKA Holder, Ø 11.5 mm	For MEDENTIKA Holder, Ø 15.8 mm
062.4601	062.4602

Variobase® XC for Crown, AH 7.0 mm

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
GH 1.5 mm	062.5017
	062.5020
GH 2.5 mm	062.5018
	062.5021
GH 3.5 mm	062.5019
	062.5022

Replacement screw

065.0036

Burn-out Copings for Variobase® XC for Crown

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
065.0106	065.0107

Variobase® XC for Crown AS, AH 7.0 mm

Ø 4.5 mm (RB/WB)	
GH 0.75 mm	—
GH 1.5 mm	062.5009
GH 2.5 mm	062.5010
GH 3.5 mm	062.5011

Replacement screw

065.0037

Burn-out Copings for Variobase® XC for Crown AS

Ø 4.5 mm (RB/WB)
065.0131

Variobase® XC for Bridge/Bar, AH 7.0 mm

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
GH 0.75 mm	—
GH 1.5 mm	062.5065
	062.5068
GH 2.5 mm	062.5066
	062.5069
GH 3.5 mm	062.5067
	062.5070

Replacement screw

065.0036

Burn-out Copings for Variobase® XC for Bridge/Bar

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
065.0131	065.0132

Variobase® XC for Bridge/Bar AS, AH 7.0 mm

Ø 4.5 mm (RB/WB)	
GH 0.75 mm	—
GH 1.5 mm	062.5073
GH 2.5 mm	062.5074
GH 3.5 mm	062.5075

Replacement screw

065.0037

Burn-out Copings for Variobase® XC for Bridge/Bar AS

Ø 4.5 mm (RB/WB)
065.0135

Variobase® C

Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
GH 1.5 mm	062.4981
	062.4982
GH 2.5 mm	062.5028
	062.5030
GH 3.5 mm	062.5029
	062.5031

Replacement screw

065.0036

Variobase® XC for Crown, AH 7 mm

Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
GH 0.75 mm	062.5023
	062.5024
GH 1.5 mm	062.5024
	062.5061
GH 2.5 mm	062.5059
	062.5062

Burn-out Copings for Variobase® XC for Crown

Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
065.0108	065.0124

Variobase® XC for Crown AS, AH 7 mm

Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
GH 0.75 mm	062.5012
	062.5013
GH 1.5 mm	062.5013
	062.5063
GH 2.5 mm	062.5060
	062.5064

Burn-out Copings for Variobase® XC for Crown AS

Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
065.0113	065.0126

Variobase® XC for Bridge/Bar, AH 7 mm

Ø 5.5 mm (WB)	
GH 0.75 mm	062.5071
	062.5072
GH 1.5 mm	062.5072

Burn-out Copings for Variobase® XC for Bridge/Bar

Ø 5.5 mm (WB)
065.0108

Variobase® XC for Bridge/Bar AS, AH 7 mm

Ø 5.5 mm (WB)	
GH 0.75 mm	062.5076
	062.5077
GH 1.5 mm	062.5077

Burn-out Copings for Variobase® XC for Bridge/Bar

Ø 5.5 mm (WB)
065.0108

Variobase® C

Ø 5.5 mm (WB)	
GH 0.75 mm	062.5032
	062.4983
GH 1.5 mm	062.4983

Replacement screw

065.0036

Novaloc®

Angle	0°	15°
1.5 mm	062.4501	—
2.5 mm	062.4502	062.4507
3.5 mm	062.4503	062.4508
4.5 mm	062.4504	062.4509
5.5 mm	062.4505	062.4510
6.5 mm	062.4506	062.4511
7.5 mm	—	062.4512

Replacement screw

065.0036

Processing Package titanium Novaloc®

2010.601-NOV

Processing Package PEEK Novaloc®

2010.611-NOV

Novaloc® Processing Insert Removal Instrument

2010.731-NOV

Novaloc® Retention Insert Instrument

2010.741-NOV

Novaloc® Matrix Housing Extraction Instrument

2010.751-NOV

Anatomic Abutments

Angle	0°	17°
2.5 mm	062.4103	062.4153
3.5 mm	062.4104	062.4154

Replacement screw

065.0036

Locator® Zest

0°	15°
GH 1 mmH	02326
GH 2 mmH	02327
GH 3 mmH	02328
GH 4 mmH	02329
GH 5 mmH	02330
GH 6 mmH	02331

Replacement screw

065.0036

Instruments and Accessories for Locator®

Up to 20° Divergence

Up to 40° Divergence

Up to 20° Divergence Locator® Zest

Up to 40° Divergence Locator® Zest

8519-10

8540-10

046.417/6

8393

Titanium Abutments
Hybrid Abutments

Angle	0°	17°
2.5 mm	062.4103	062.4153
3.5 mm	062.4104	062.4154

Replacement screw

065.0036