

Slim emergence profile (RB/WB) prosthetics

fit on all BLC and BLX implants (Ø 3.3 mm - Ø 6.5 mm)

- crowns with slim emergence profiles
- bridges

Implant-abutment compatibility		
RB/WB and WB		RB Implants WB Implants
«RB/WB» RB/WB Variobase® XC for Crown	compatible with	
«WB» WB Variobase® XC for Crown	compatible with	No WB abutments on RB implants!

Wide Emergence Profile (WB) prosthetics

fit only on WB implants (BLC Ø 4.5 mm - Ø 6.5 mm, BLX Ø 5.0 mm - Ø 6.5 mm)

- crowns with wide emergence profiles

Abutment Level, on Screw-retained Abutments, Ø 4.6 mm		
Implant Analog for Screw-retained Abutments, Ø 4.6 mm, 0°	Open-tray Impression for Screw-retained Abutments, Ø 4.6 mm	Closed-tray Impression for Screw-retained Abutments, Ø 4.6 mm
023.4756	Crown 025.2244 Bridge/Bar 025.0012	Crown 025.2246 Bridge/Bar 025.0014
Scanbody for Screw-retained Abutments, Ø 4.6 mm	Repositionable Implant Analog for Screw-retained Abutments, Ø 4.6 mm	RevEX SRA 4.6
025.0081	025.0008	025.0083

Implant Level			
Analog	Open-tray	Closed-tray	
065.0021	065.0031	065.0033	065.4310
Straumann® Scanbody	Repositionable Implant Analog		
065.0103	065.0023		

Temporary Titanium Coping for Screw-retained Abutments, Ø 4.6 mm		
Titanium Copings	Ø 4.6 mm (RB/WB)	
	Crown	023.4747
	Bridge	023.4751
incl. screw		023.4763
	Replacement screw	023.4763
	SRA Protective Caps	
	Ø 4.6 mm	
	5.1 mm	024.4323-04
	6.6 mm	024.4324-04
	8.1 mm	024.4325-04

RB/WB Temporary Abutment, VITA CAD-Temp®		
	Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
	GH	1.5 mm
		064.4361
		064.4371

Temporary Abutments for Crown, AH 10 mm		
	Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
	GH	1.5 mm
		064.4362
		064.4372
		064.4382
		064.4373
		064.4383
		065.0036

Temporary Abutments for Bridge, AH 10 mm		
	Ø 4.5 mm (RB/WB)	
	GH	1.5 mm
		064.4352
		065.0036

Temporary Abutments for Crown, AH 11 mm		
	Ø 5.5 mm (WB)	
	GH	0.75 mm
		064.4391
		064.4392
		065.0036

WB Temporary Abutment, VITA CAD-Temp®		
	Ø 5.5 mm (WB)	
	GH	1.5 mm
		064.4390

Screw-retained Abutments

Pre-milled Abutment Blanks

Variobase®

Screw-retained Abutments		
	Ø 4.6 mm (RB/WB)	
	Angle	0°
		17°
		30°
	0.75 mm	—
	1.5 mm	062.47225
	2.5 mm	062.47235
	3.5 mm	062.47245
	4.5 mm	062.47255
	5.5 mm	062.47355
		062.47435
		062.47445
		062.47455
		065.0036

RB/WB SRA PLAN Abutment, POM		
	0°	17°
		30°
		025.0073V4
		025.0074V4
		025.0075V4

RB/WB Pre-milled Abutment Blanks		
For MEDENTIKA Holder, Ø 11.5 mm	For MEDENTIKA Holder, Ø 15.8 mm	
		062.4601
		062.4602

Variobase® XC for Crown, AH 7.0 mm		
	Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
	GH	1.5 mm
		062.5017
		062.5020
		062.5018
		062.5021
		062.5019
		062.5022
		065.0036

Variobase® XC for Crown AS, AH 7.0 mm		
	Ø 4.5 mm (RB/WB)	
	GH	0.75 mm
		—
		062.5009
		062.5010
		062.5011
		065.0037

Variobase® XC for Bridge/Bar, AH 7.0 mm		
	Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
	GH	0.75 mm
		—
		062.5065
		062.5068
		062.5066
		062.5069
		062.5067
		062.5070
		065.0036

Variobase® XC for Bridge/Bar AS, AH 7.0 mm		
	Ø 4.5 mm (RB/WB)	
	GH	0.75 mm
		—
		062.5073
		062.5074
		062.5075
		065.0037

Variobase® C		
	Ø 3.8 mm (RB/WB)	Ø 4.5 mm (RB/WB)
	GH	1.5 mm
		062.4981
		062.4982
		062.5028
		062.5030
		062.5029
		062.5031
		065.0036

Variobase® XC for Crown, AH 7 mm		
	Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
	GH	0.75 mm
		062.5023
		062.5024
		062.5061
		062.5059
		062.5062

Burn-out Copings for Variobase® XC for Crown		
	Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
		065.0108
		065.0124

Variobase® XC for Crown AS, AH 7 mm		
	Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
	GH	0.75 mm
		062.5012
		062.5013
		062.5063
		062.5060
		062.5064

Burn-out Copings for Variobase® XC for Crown AS		
	Ø 5.5 mm (WB)	Ø 6.5 mm (WB)
		065.0113
		065.0126

Variobase® XC for Bridge/Bar, AH 7 mm		
	Ø 5.5 mm (WB)	
	GH	0.75 mm
		062.5071
		062.5072

Burn-out Copings for Variobase® XC for Bridge/Bar		
	Ø 5.5 mm (WB)	
		065.0108

Variobase® XC for Bridge/Bar AS, AH 7 mm		
	Ø 5.5 mm (WB)	
	GH	0.75 mm
		062.5076
		062.5077

Burn-out Copings for Variobase® XC for Bridge/Bar		
	Ø 5.5 mm (WB)	
		065.0108

Variobase® C		
	Ø 5.5 mm (WB)	
	GH	0.75 mm
		062.5032
		062.4983
		065.0036

Anatomic Abutments		
	Ø 6.0 mm (RB/WB)	
	Angle	0°
		17°
	2.5 mm	062.4103
		062.4153
	3.5 mm	062.4104
		062.4154
		065.0036

Locator® Zest		
	0°	15°
	GH	
	1 mmH	02326
	2 mmH	02327
	3 mmH	02328
	4 mmH	02329
	5 mmH	02330
	6 mmH	02331
		02398A-SB
		02399A-SB
		02400A-SB
		02401A-SB
		02402A-SB
		02403A-SB

Instruments and Accessories for Locator®		
Up to 20° Divergence	Up to 40° Divergence	
		8519-10
		8540-10
		046.417/6
		8393

Up to 20° Divergence Locator® Zest		
		8519-10

Up to 40° Divergence Locator® Zest		
		8540-10

Novaloc®		
	Ø 3.8 mm (RB/WB)	
	Angle	0°
		15°
	1.5 mm	062.4501
	2.5 mm	062.4502
	3.5 mm	062.4503
	4.5 mm	062.4504
	5.5 mm	062.4505
	6.5 mm	062.4506
	7.5 mm	062.4512
		065.0036

Processing Package titanium Novaloc®		
		2010.601-NOV

Processing Package PEEK Novaloc®		
		2010.611-NOV

Novaloc® Processing Insert Removal Instrument		
		2010.731-NOV

Novaloc® Retention Insert Instrument		
		2010.741-NOV

Novaloc® Matrix Housing Extraction Instrument		
		2010.751-NOV